



PHOTOCATHODE MOUNTING ON MOLYBDENUM PUCK BODY

CIS

Jefferson Lab

This note summarizes the steps to mount photocathode to molybdenum puck body described in drawing JL0047248.

1. Collect items
 - a. Mo puck, photocathode removed
 - b. Tantalum cup, choose ID
 - c. Photocathode (nominal 600um thick)
 - d. Ta cup recessed pressing tool
 - e. Ta cup crimping tool
 - f. Indium foil
 - g. Indium tweezers (only touch indium foil)
 - h. Non-indium tweezers (only touch photocathode or moly puck)
 - i. Pliers to lift puck from hot plate
2. Dry fit tests
 - a. Place photocathode on puck recess, ensure it fits, then remove
 - b. Place Ta cup on puck, ensure fit is not too tight or loose, carefully use Ta cup pressing tool if helpful, but do not dent Ta cup
 - c. Carefully remove Ta cup
3. Mount photocathode under GN2 glove box
 - a. Cut indium square, place on puck center
 - b. Place puck on heater, start heater, indium will melt 156 C. Note - thermocouples typically read over 200C due to temperature gradients
 - c. Place photocathode on indium, press flat on edges
 - d. Check that puck is well seated, flat; excess indium is OK, remove later
 - e. Remove puck from heater using pliers, to cooling plate
4. Mount Ta cup (once cool, outside glove box)
 - a. If indium 'ball' remove w/ vertical force, don't push photocathode
 - b. Place Ta cup on puck
 - c. Press down with recessed pressing tool if required/desired
 - d. Place crimping screws on puck, thumb screw until uniform circle
 - e. Place crimping cylinder over puck/screw
 - f. Test Ta cup well seated (add/remove crimping cylinder)
 - g. Turn opposing pairs of thumb screws 1 turn at a time, until Ta cup is sufficiently crimped to puck body; note – turning until stopped will make Ta cup removal difficult
 - h. Remove Ta crimping tool
 - i. Use Ta cup recessed pressing tool if desired