



Ultra Detectors

Storage, Handling, and Cleaning Instructions

- Store the detectors at room temperature in clean air.
- Protect them from dust and corrosive fumes. The corrosive fumes will damage the external as well as the internal parts. Electrical contacts are very critical to the normal operation of the detector. Corrosion can form high-resistance paths in these contacts and ruin the detector.
- Avoid mechanical shocks and protect from mechanical damage.
- Wear protective latex or polyethylene gloves when handling the detectors. Finger grease and body salts can destroy the integrity of electrical contacts.
- Keep the protective cap in place when the detector is not in use. This will protect the fragile front surface from accidental damage. This is especially important when the detector is being installed, removed, or transported to other sites.
- After prolonged use, it is normal for the detector to require cleaning of the front, critical surface, especially if the detector surface is suspected of being contaminated, causing higher background. Use “fluffed” cotton swabs wetted with methanol to gently rub the surface. Change the swabs frequently to avoid redepositing the contaminants. **DO NOT FLOOD THE DETECTOR SURFACE WITH METHANOL. THE CONDUCTIVE ELASTOMER CAN SWELL AND ELECTRICAL CONTACT COULD BE LOST, RESULTING IN NOISE AND SEVERELY COMPROMISED RESOLUTION.**
- Blow clean, dry nitrogen to dry the surface. Do not use filter paper to dry the surface. Although the implanted contact is “buried” under the surface, it is very shallow and can be easily damaged. **The slightest damage to even a small area will ruin the detector.**

If a detector must be sent back to ORTEC, it should be properly packaged with the PROTECTIVE CAP installed. Use the packaging that the detector was shipped in, to avoid damage during transit. Improper handling or packaging may void the warranty.

804558B / 0405
Printed in USA



Advanced Measurement Technology, Inc.

ORTEC®


801 South Illinois Avenue
Oak Ridge, TN 37831-0895
Tel. (865) 482-4411 • www.ortec-online.com

AMETEK®

Quality Assurance Data Sheet

ULTRA™ Alpha Detector

Performance Specifications

Model No.	BU-013-050-1000-S	Serial No.	40-015 A12
Date Shipped	6-13-07		
Bias Voltage	115	volts (positive)	
Noise	6.9	keV (FWHM)	
Warranted Alpha Resolution	13	keV (FWHM)	
Special Notes (if any)			
Data Certified By:  [signature]			

Handling

Wear protective latex or polyethylene gloves when handling the detectors. Finger grease and body salts can destroy the integrity of the electrical contacts.

The ULTRA series detector has a thin (500 Å), ion-implanted contact immediately under the silicon surface. The contact is, thus, more rugged than that formed by an evaporated gold layer. If the silicon surface is scratched, however, the detector will be damaged. Reasonable precautions should therefore be taken when handling these detectors.

If a detector must be sent back to ORTEC, it should be properly packaged with the **PROTECTIVE CAP** installed. Use the packaging that the detector was shipped in to avoid damage during transit. Improper handling or packaging may void the warranty.

Application of Bias

The detector connector's central electrode should be provided with positive bias and produces a negative output signal. The recommended operating voltage is given in the performance specifications. This bias voltage is carefully selected to provide optimum field strength and should not be exceeded.

Environmental Considerations

Although there are no epoxies or rings to absorb moisture, it is essential that the detector be dry during operation. Moisture can cause both high leakage current and low breakdown voltage. To remove condensed water or condensed organic vapors, the detector should be placed in vacuum for a short time before applying bias.

Cleaning

If contaminated, the detector surface can be cleaned with a cotton swab moistened with methanol. Gently rub the detector surface with the swab. Repeat with fresh methanol. Blow dry with dry nitrogen gas. Before applying bias, leave the detector under vacuum for 30 minutes to remove all moisture from the surface.

Measurement of Alpha Resolution

The measurement of alpha particle resolution must be performed in vacuum with a uniform, ultra-thin source located at a source-to-detector distance of 1-1/2 times the detector diameter.



Advanced Measurement Technology, Inc.

ORTEC®


801 South Illinois Avenue
Oak Ridge, TN 37831-0895
Tel. (865) 482-4411 • www.ortec-online.com



Quality Assurance Data Sheet

ULTRA™ Alpha Detector

Performance Specifications

Model No.	BU-013-050-1000-S	Serial No.	40-015 A11
Date Shipped	6-13-07		
Bias Voltage	115	volts (positive)	
Noise	6.9	keV (FWHM)	
Warranted Alpha Resolution	13	keV (FWHM)	
Special Notes (if any)			
Data Certified By:  [signature]			

Handling

Wear protective latex or polyethylene gloves when handling the detectors. Finger grease and body salts can destroy the integrity of the electrical contacts.

The ULTRA series detector has a thin (500 Å), ion-implanted contact immediately under the silicon surface. The contact is, thus, more rugged than that formed by an evaporated gold layer. If the silicon surface is scratched, however, the detector will be damaged. Reasonable precautions should therefore be taken when handling these detectors.

If a detector must be sent back to ORTEC, it should be properly packaged with the **PROTECTIVE CAP** installed. Use the packaging that the detector was shipped in to avoid damage during transit. Improper handling or packaging may void the warranty.

Application of Bias

The detector connector's central electrode should be provided with positive bias and produces a negative output signal. The recommended operating voltage is given in the performance specifications. This bias voltage is carefully selected to provide optimum field strength and should not be exceeded.

Environmental Considerations

Although there are no epoxies or rings to absorb moisture, it is essential that the detector be dry during operation. Moisture can cause both high leakage current and low breakdown voltage. To remove condensed water or condensed organic vapors, the detector should be placed in vacuum for a short time before applying bias.

Cleaning

If contaminated, the detector surface can be cleaned with a cotton swab moistened with methanol. Gently rub the detector surface with the swab. Repeat with fresh methanol. Blow dry with dry nitrogen gas. Before applying bias, leave the detector under vacuum for 30 minutes to remove all moisture from the surface.

Measurement of Alpha Resolution

The measurement of alpha particle resolution must be performed in vacuum with a uniform, ultra-thin source located at a source-to-detector distance of 1-1/2 times the detector diameter.

RTEC®
Princeton Applied Research
SIGNAL RECOVERY

Clartron Analytical

E.I.N. #62-1369733

**Advanced
Measurement
Technology, Inc.**

801 South Illinois Avenue
Oak Ridge, TN 37831-0895
Tel: (865) 482-4411
Fax: (865) 483-2279

SALES ORDER NO.	PRINTED ON	PAGE
0171459	06/11/07	
ORDER DATE	EXPORT LICENSE OR LETTER OF CREDIT	
05/16/07		
CUSTOMER PURCHASER ORDER NUMBER		
07-P1139		

SHIPPING DOCUMENT

SHIPMENT NO.	DATE SHIPPED	BILL OF LADING NO.
0174002		271084
METHOD OF SHIPMENT	FREIGHT TERMS	WAREHOUSE
BEST WAY	FOB DEST	USM

JEFFERSON LAB

12000 JEFFERSON AVE.

NEWPORT NEWS, VA 23606

BILL
TO

JEFFERSON SCIENCE ASSOCIATES
(JSA) JEFFERSON LAB (JLAB)
ATTN: ACCT. PAYABLE
628 HOFSTADTER ROAD
NEWPORT NEWS, VA
23606

SHIP TO CUSTOMER

1TJNAF-01

BILL TO CUSTOMER

1TJNAF

THESE COMMODITIES ARE LICENSED BY THE UNITED STATES GOVERNMENT FOR THE COUNTRY OF ULTIMATE DESTINATION THAT IS INDICATED IN THE "SHIP TO ADDRESS" UNLESS OTHERWISE NOTED. RE-EXPORT OF THESE PRODUCTS REQUIRES PRIOR WRITTEN AUTHORIZATION FROM THE UNITED STATES GOVERNMENT. DIVERSION OR UNAUTHORIZED RE-EXPORT OF THESE COMMODITIES IS CONTRARY TO UNITED STATES LAW AND IS STRICTLY PROHIBITED.

LINE	PRODUCT NUMBER	PRODUCT DESCRIPTION/SERIAL NUMBER	CURRENT QTY. BACKORDERED	QTY. SHIPPED	NEW QTY. BACKORDERED
		PRICE INCLUDES 1.5% FOR DESTINATION FEE.			
1	142A	PREAMPLIFIER	2.00	2	
2	U-OFSTD	7152253 715225 OFFSTANDARD ULTRA DETECTOR	2.00	2.0	4
		BU-013-050-1000-S			
		40-015A11, 40-015A12			

CUSTOMER PACKING LIST

