

Equipment

- Two 1" diameter (or greater) $f=1\text{m}$ lenses, each with a mount which controls both angle and lens position (i.e. a 4axis or 5 axis lens mount) + stands for them
- The quad photodiode
- Control over the helicity board
- 3-5mW of Hall A laser for alignment (CW or pulsed, either is fine)
- $>20\mu\text{A}$ of electron beam (preferably $\sim 70\mu\text{A}$, but not strictly necessary)
- (Conditions of injector beamline should be as if accelerator were going to run $70\mu\text{A}$ of 1GeV beam, or $150\mu\text{A}$ of 2GeV beam)

People – Caryn Palatchi, Sachinthani Premathilake, Ciprian Gal, Kent Paschke

Need walk-through of injector laser for Sachinthani Premathilake

- Day1 benchmarking
 - **Morning**
 - *Need someone in control room who can change beam current, turn on/off autogaining on bpms*
 - HallA Electron beam $>20\mu\text{A}$ ($70\mu\text{A}$ is good) going up to at least FC1
 - BCM/BPM calibration scan - $5\mu\text{A}$ steps of current up to max current, auto gaining on injector bpms off
 - Autogaining of injector bpms back on
 - **Afternoon**
 - 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
 - Tweak Pockels cell translation – 1 hour
 - Repeat 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
 - **Evening**
 - *Need Access to injector laser room*
 - *Need someone who can get the laser to give us 3-5mW of Hall A laser beam*
 - Get spiricon measure of spot size at cathode
 - Repeat measure of spot size at pockels cell (will bring our own spiricon for this)
-
- Day2 downstream lens insertion
 - **Morning**
 - *Need Access to injector laser room*
 - *Need someone who can get the laser to give us 3-5mW of Hall A laser beam*
 - 3-5mW Hall A laser (CW or pulsed, either is fine)
 - Insert 1m lens downstream of Pockels Cell at predetermined z-position (personal preference is upstream of clean-up polarizer, but only Hall A laser spot size will be affected)
 - Measure spiricon spot size at cathode
 - **Afternoon**

- *Need someone in control room who can change beam current*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
- **Evening**
- *Need Access to injector laser room*
- *Need someone who can get the laser to give us 3-5mW of Hall A laser beam*
- 3-5mW Hall A laser (CW or pulsed, either is fine)
- Start Setup of QPD pickoff / calibration

- Day 3 upstream lens insertion and PC re-alignment
- **Morning**
- *Need Access to injector laser room*
- *Need someone who can get the laser to give us 3-5mW of Hall A laser beam*
- 3-5mW Hall A laser (CW or pulsed, either is fine)
- Remove downstream 1m lens
- Insert 1m lens upstream of Pockels Cell at predetermined z-position
- measure of spot size at pockels cell (will bring our own spiricon for this)
- Measure divergence of laser at Pockels cell
- Measure spiricon spot size at cathode
- Setup pick off to QPD
- Measure spot size at QPD
- **Afternoon**
- Calibrate QPD
- Check PC alignment starting point– S1, S2, no anal, RHWP scan
- Align Pockels cell
- **Evening**
- *Need someone in control room who can change beam current*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
- PC translation if needed
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours

- Day 4 PC-realignment (if needed) + downstream lens insertion
- **Morning**
- *Need Access to injector laser room*
- *Need someone who can get the laser to give us 3-5mW of Hall A laser beam*
- 3-5mW Hall A laser (CW or pulsed, either is fine)
- Insert 1m lens downstream of Pockels Cell at predetermined z-position
- Measure spiricon spot size at cathode

- **Afternoon**
- *Need someone in control room who can change beam current, turn on*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
- *Need Access to injector laser room*
- *Need someone who can get the laser to give us 3-5mW of Hall A laser beam*
- Remove 1m lens downstream of PC
- **Evening**
- Allow time opportunity for any PC realignments that seem necessary on the QPD
- *Need someone in control room who can change beam current*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours

- Day 5 Photocathode rotation and PC translation
- **Morning**
- *Need Access to injector room*
- *Need someone who can help us rotate the photocathode Angle #2*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
- **Afternoon**
- *Need Access to injector room*
- *Need someone who can help us rotate the photocathode Angle #3*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
- *Need Access to injector room*
- *Need someone who can help us rotate the photocathode Angle #4*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
- **Evening**
- *Need Access to injector room*
- *Need someone who can help us rotate the photocathode FINAL ANGLE*
- HallA Electron beam >20uA (70uA is good) going up to at least FC1
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours
- PC translation to optimize
- 4 RHWP scans (IHWP in/out PITA 0/ PITA non-zero) – 2 hours